

SELF-EVALUATION:
 1=Independent
 2= Need assist

Name: _____ ID: _____ Department: _____

Competency: **Staff will properly use the Accucheck Inform® to obtain Point-of Care Glucose measurement**

VALIDATOR EVALUATION:
 1= Independent W/O Errors
 2= Independent w/ Minor Errors
 3= Verbal Prompts w/ Errors

Skill/Steps	Rationale	Self		Validator				
		1	2	1	2	3		
Performs high and low quality control tests, following Infection Control guidelines 1. Checks expirations dates of supplies, discarding outdated 2. Performs QC each day or when meter error is suspected 3. Turns meter on and scans own OPERATOR ID barcode (i.e. UMID#). 4. Selects CONTROL TEST then Level 1 or 2 5. Scans & verifies lot number of Control Solutions 6. Scans test strip vial code & verifies match on screen 7. Removes strip, recaps vial immediately, inserts strip into when indicated on display 8. Applies 1 drop control solution to strip's yellow dosing area when indicated 9. Notes result, runs next control level a. FAIL: enters appropriate comment(s), correct potential causes, repeats control test b. Continued "out of range", reports to Charge Nurse, removes from service, exchanges meter at Chemistry Dept. or local lab	Annual re-certification may be done and recorded electronically through glucose meter 1. Control solutions are to be relabeled with current date & 3 month expiration date upon opening. Use until that or printed expiration date, whichever occurs first. DO NOT SHAKE CONTROL SOLUTIONS Strips expire on manufacturer's date if tightly closed in original vial 2. QCs prior to testing patient samples 3., 5, 6. Manual entry of data increases risk of error Point of Care Policy Accu-Chek Inform II Quick Reference Guide AccuChek Inform II Operator's Manual	1	2	1	2	3		
Comments:								
Accurately measures and records capillary blood glucose reading following Infection Control guidelines 1. Checks expirations dates of supplies, discarding outdated 2. Identifies patient with 2 identifiers 3. Instructs patient re: procedure and to wash hands, if able 4. Turns meter on and scans own OPERATOR ID barcode (i.e. UMID #) 5. Ensures wireless connection is accessed, if available 6. Selects PATIENT TEST 7. Scans patient's bar coded wristband at bedside & verifies that PATIENT ID# on screen matches wristband's CSN 8. Scans then verifies that test strip code on vial and meter match 9. Removes strip, recaps vial immediately, inserts strip into meter when displayed on screen 10. When indicated, obtains pure blood capillary, venous, or arterial sample: a. Allows alcohol to dry, if used for prep b. Gently squeezes, pricks side of fingertip or heel with lancet c. Wipes away first drop of blood with dry gauze	This may be completed with a test reagent or patient sample. Current Accu-chek Inform users may describe steps 1. Strips expire on manufacturer's date if tightly closed in original vial 4., 7, 8. Manual entry of data increases risk of error. 7. Carefully enters CSN manually, only in areas that do not use ID band MRN is not correct entry for glucose meter 10. Prevents error readings Note- CORD Blood is not acceptable	Self	Validator	1	2	1	2	3
Comments:								

Competency

Employee Name

UMID

<p>Accurately measures and records capillary blood glucose reading following Infection Control guidelines</p>	<p><i>This may be completed with a test reagent or patient sample. Current Accu-chek Inform users</i></p>	<p>Self</p>		<p>Validator</p>		
<ol style="list-style-type: none"> 1. Checks expirations dates of supplies, discarding outdated 2. Identifies patient with 2 identifiers 3. Instructs patient re: procedure and to wash hands, if able 4. Turns meter on and scans own OPERATOR ID barcode (i.e. UMID #) 5. Ensures wireless connection is accessed, if available 6. Selects PATIENT TEST 7. Scans patient's bar coded wristband at bedside & verifies that PATIENT ID# on screen matches wristband's CSN 8. Scans then verifies that test strip code on vial and meter match 		<p>1</p>	<p>2</p>	<p>1</p>	<p>2</p>	<p>3</p>
		<p>Comments:</p>				

This employee has satisfactorily completed the above competency

<p>Employee</p>	<p>Signature</p>	<p>initials</p>	<p>Date</p>
<p>Validator</p>	<p>Signature</p>	<p>Initials</p>	<p>Date</p>

This competency has not been achieved. Action Plan and expected date of completion as follows:

<p>Validator</p>	<p>Signature</p>	<p>Initials</p>	<p>Date</p>
------------------	------------------	-----------------	-------------